

Optional accessories for autoclave

Wire basket


The structure easily intaking the steam.
Selectable depending on the diameter and the height of chambers.



Applicable model	Product code	Inner size (Dia.xHeight)
MC-23, TR-24S, TR-24LA/B, KTR-2322, KTR-2346A/B	BW2120	φ210×200mm
MC-23, TR-24LA/B, KTR-2346A/B	BW2140	φ210×400mm
KTR-3022, KTR-3045A/B, KTR-3065A/B, MC-30S, MC-30L, MC-3032S, MC-3032L, RK-3030	BW2720	φ270×200mm
KTR-3045A/B, KTR-3065A/B, MC-30S, MC-30L, MC-3032S, MC-3032L, RK-3030	BW2730	φ270×300mm
KTR-3045A/B, KTR-3065A/B, MS-30S, MC-30L, MC-3032S, MC-3032L, RK-3030	BW2740	φ270×400mm
CLG-32L, CLG-32LDVP, CL-32S(DP), CL-32L(DP), CLM-32L	BW3020	φ300×200mm
CLG-32L, CLG-32LDVP, CL-32S(DP), CL-32L(DP), CLM-32L	BW3030	φ300×300mm
CLG-40M(DVP), CLG-40L(DVP), CLS-40S(DP), CLS-40M(DP), CLS-40L(DP), CLM-40M, CLM-40L, RKZ-40II	BW3815	φ380×150mm
CLG-40M(DVP), CLG-40L(DVP), CL-40S(DP), CL-40M(DP), CL-40L(DP), CLM-40M, CLM-40L, RKZ-40II	BW3820	φ380×200mm
CLG-40M(DVP), CLG-40L(DVP), CL-40S(DP), CL-40M(DP), CL-40L(DP), CLM-40M, CLM-40L, RKZ-40II	BW3830	φ380×300mm
CLG-40M(DVP), CLG-40L(DVP), CLM-40M, CLM-40L, CL-40S(DP), CL-40M(DP), CL-40L(DP)	BW3840	φ380×400mm

Bottom Pan for wire basket

Enables to prevent spills of liquid culture media during sterilization. As your baskets can store these pans, you can also purchase it seperately as necessary.



Applicable wire basket (Product code)	Product code	Inner size (Dia.xHeight)
BW2120, BW2140	DH2036	φ203×60mm
BW2720, BW2730, BW2740	DH2676	φ267×60mm
BW3020, BW3030	DH2896	φ289×60mm
BW3820, BW3830, BW3840	DH3636	φ363×60mm

Bucket

It is useful for sterilizing a large quantity of soil and liquid that wire basket is not suitable for. Also, it is recommended to use for sterilizing a waste when not using a sterilization bag.



Applicable model	Product code	Inner size (Dia.xHeight)
MC-23, TR-24LA/B, KTR-2346A/B	VN2121	φ210×210mm
MC-30S, MC-30L, KTR-3045A/B, KTR-3065A/B, CLG-32L(DVP), CL-32S(DP), CL-32L(DP), CLM-32L	VN2727	φ270×270mm
CLG-40M(DVP), CLG-40L(DVP), CLM-40M, CLM-40L, CL-40S(DP), CL-40M(DP), CL-40L(DP)	VN3630	φ363×300mm

Do not stack the buckets in two layers (upper & lower) in the chamber. Otherwise, the air in the lower bucket is not evacuated, which causes a sterilization failure.

Dressing drum

Suitable to store clothes. The slits at open makes the objects exposed with steam and that at close prevents bacteria invaded after sterilization.



Applicable model	Product code	Inner size (Dia.xHeight)
All Models	SC2121	φ180×120mm
All Models except for MC-23, TR-24S, TR-24LA/B, KTR-2322, KTR-2346A/B	SC2727	φ270×180mm

The optional accessories above are available both when you place an order for the machine and after the order.