

## Información general

**Nombre de producto:** ASPIRADORES DE PIPETAS DE SEGURIDAD

**Fecha Rev.:** Jul 11, 2023

**Detalles:**

Acepta pipetas estándar y desechables.

Evita riesgos incluso con líquidos tóxicos o corrosivos.

Debido a su alojamiento flexible, las pipetas se ajustan suavemente.

Fácil manejo con una sola mano. Girando la rueda dentada hacia arriba o hacia abajo, se obtiene un llenado o vaciado preciso. Pulsando la clavija lateral se produce un vaciado automático. Fácilmente desmontable.

## Fotografía



## Datos técnicos

## Ficha de producto

Ref.	Rango (mL)	Color
10ASPS0102	0-2	Azul
10ASPS0110	0-10	Verde
10ASPS0125	0-25	Rojo



## General Information

**Product name:** PIPPETE PUMPS

**Rev. Date:** Jul 11, 2023

**Details:**

For standard and disposable pipettes.

Avoid risks even with corrosive and toxic liquids.

Due to its flexible hole, pipettes fit in smoothly.

Convenient single-handed operation. With forward and backward wheel movements, a precise filling or emptying can be achieved. With a lever for quick liquid release. Easy disassembly.

## Picture



## Technical Information

## Product datasheet

Ref.	Range (mL)	Colour
10ASPS0102	0-2	Blue
10ASPS0110	0-10	Verde
10ASPS0125	0-25	Red

