

## Información general

**Nombre de producto:** SCHILLING BURETTES, CLASS A, COMPLETE SET - AMBER

**Fecha Rev.:** Jul 11, 2023

**Detalles:**

- Con enrase automático a cero.
- Circuito cerrado para protección de la contaminación exterior.
- Unidad compacta que hace que no sean necesarios accesorios metálicos (pinzas, etc.).
- Fácil carga de la bureta, simplemente por presión del frasco.

## Fotografía



## Datos técnicos

AMBER

## Ficha de producto

Ref.	Llave	Capacidad (ML)	Tolerancia ( $\pm$ ML)	Subdivisiones (ML)
0100491478AM	PTFE Aguja	25	$\pm 0.05$	0.1
0100491479AM	PTFE Aguja	50	$\pm 0.05$	0.1
0100491481AM	Plástico	25	$\pm 0.05$	0.1
0100491482AM	Plástico	50	$\pm 0.05$	0.1



## General Information

**Product name:** SCHILLING BURETTES, CLASS A, COMPLETE SET - AMBER

**Rev. Date:** Jul 11, 2023

**Details:**

- With automatic zeroing.
- Closed circuit to avoid external contamination.
- Compact unit with no need for metallic accesories (clamps, etc.).
- Easy loading from flask pressure.

## Picture



## Technical Information

AMBER

## Product datasheet

Ref.	Stopcock	Capacity (ML)	Tolerance ( $\pm$ ML)	Subdivisiones (ML)
0100491478AM	PTFE NEEDLE	25 ML	$\pm 0.05$	0.1
0100491479AM	PTFE NEEDLE	50 ML	$\pm 0.05$	0.1
0100491481AM	Plastic	25 ML	$\pm 0.05$	0.1
0100491482AM	Plastic	50 ML	$\pm 0.05$	0.1

