

## Información general

**Nombre de producto:** DESECADORES EN VIDRIO

**Fecha Rev.:** Jul 11, 2023

**Detalles:**

Vidrio borosilicato, tapa botón a rosca, sin placa.

## Fotografía



## Datos técnicos

| Ref.       | Ø Interior (MM) | Capacidad Aprox. (ML) |
|------------|-----------------|-----------------------|
| 0200929210 | 150             | 1500                  |
| 0200929220 | 200             | 4000                  |
| 0200929250 | 250             | 6000                  |
| 0200929260 | 300             | 10000                 |

## General Information

**Product name:** GLASS DESICCATORS

**Rev. Date:** Jul 11, 2023

**Details:**

Borosilicate glass, button screw cover, without plate.

## Picture



## Technical Information

| Ref.       | Ø ID (MM) | Capacity Approx. (ML) |
|------------|-----------|-----------------------|
| 0200929210 | 150       | 1500                  |
| 0200929220 | 200       | 4000                  |
| 0200929250 | 250       | 6000                  |
| 0200929260 | 300       | 10000                 |