

Información general

Nombre de producto: EQUIPO DE FILTRACIÓN CON ERLLENMEYER Y SOPORTE FILTRANTE CON OLIVA, UNIÓN ESMERILADA

Fecha Rev.: Jul 11, 2023

Detalles:

Compuesto por los siguientes productos:

- Matraz Erlenmeyer de 1000mL (Ref. 0100259910) y 2000mL (Ref. 0100259912)
- Cuerpo con placa porosa, hembra. Con oliva de vidrio (Ref. 0100259933) y plástico (Ref. 0100259934)
- Embudo de filtración de 300mL (Ref. 0100260004)
- Pinza de sujeción (Ref. 0100260006)

Fotografía



Datos técnicos

Ficha de producto

Ref.	Oliva	Capacidad (ML)
0100259902	Plástico	2000
0100259903	Vidrio	1000
0100259904	Plástico	1000



General Information

Product name: FILTRATION SET WITH ERLLENMEYER AND FILTER HOLDER WITH OLIVE, GROUND JOINT

Rev. Date: Jul 11, 2023

Details:

Includes the following products:

- Erlenmeyer flask. 1000mL (Ref. 0100259910) and 2000mL (Ref. 0100259912)
- Body with fritted disc, socket joint. With glass (Ref. 0100259933) and plastic olive (Ref. 0100259934)
- Filter funnel, 300mL (Ref. 0100260004)
- Clamp (Ref. 0100260006)

Picture



Technical Information

Product datasheet

Ref.	Hose	Capacity (ML)
0100259902	Plastic	2000 ML
0100259903	Glass	1000 ML
0100259904	Plastic	1000 ML

The logo consists of the word "LORAMEL" written vertically in a large, light green, sans-serif font. The letters are bold and have a slight shadow effect, giving them a three-dimensional appearance. The 'L' is at the top, followed by 'O', 'R', 'A', 'M', 'E', and 'L' at the bottom.