

Información general

Nombre de producto: MATRACES AFORADOS DE KOHLRAUSCH CLASE A

Fecha Rev.: Jul 11, 2023

Detalles:

Con cuello ensanchado para una fácil introducción de líquidos viscosos y muestras sólidas. Tamaños de 100 y 200 mL típicamente usados para análisis de azúcar y de viscosidad de aceite por el método de Engler, respectivamente.

Fabricado en vidrio borosilicato 3.3.

Fotografía



Datos técnicos

Ref.	Capacidad (ML)	Tolerancia (\pm ML)
0100496910	100	± 0.100

Ficha de producto

Ref.	Capacidad (ML)	Tolerancia (\pm ML)
0100496920	200	± 0.150

MEMO

General Information

Product name: VOLUMETRIC FLASKS, ACC. TO KOHLRAUSCH, CLASS A

Rev. Date: Jul 11, 2023

Details:

With an expanded neck for an easy introduction of viscous liquids or solid samples. 100 and 200 mL sizes typically used for sugar analysis and determination of oil viscosity by the Engler method, respectively.

Made from borosilicate glass 3.3.

Picture



Technical Information

Ref.	Capacity (ML)	Tolerance (\pm ML)
0100496910	100 ML	0.10 ML

Product datasheet

Ref.	Capacity (ML)	Tolerance (\pm ML)
0100496920	200	0.15 ML

