

Información general

Nombre de producto: MATRACES PARA CULTIVO ERLLENMEYER CON 3 HENDIDURAS

Fecha Rev.: Jul 11, 2023

Detalles:

- Fabricado a partir de vidrio borosilicato 3.3.
- Presenta 3 deflectores para aportar agitación óptima y constante.

Fotografía



Datos técnicos

Ref.	Capacidad (ML)
0100126510	100
0100126525	250

Ficha de producto

Ref.	Capacidad (ML)
0100126550	500
0100126570	1000
0100126580	2000



General Information

Product name: ERLLENMEYER FLASKS FOR CULTURE WITH 3 BAFFLES

Rev. Date: Jul 11, 2023

Details:

- Manufactured from borosilicate glass 3.3.
- Features 3 baffles for optimal and constant agitation.

Picture



Technical Information

Ref.	Capacity (ML)
0100126510	100 ML
0100126525	250 ML

Product datasheet

Ref.	Capacity (ML)
0100126550	500
0100126570	1000 ML
0100126580	2000 ML

A large, light green watermark logo is oriented vertically on the right side of the page. The logo consists of the letters 'L', 'O', 'R', 'A', 'M', 'E', 'L' stacked from top to bottom. The letters are bold and stylized, with the 'L' at the top, followed by 'O', 'R', 'A', 'M', 'E', and 'L' at the bottom. The 'L' is a simple vertical bar with a horizontal top bar. The 'O' is a circle with a horizontal line through its center. The 'R' is a vertical bar with a curved top and a horizontal base. The 'A' is a vertical bar with a triangular top. The 'M' is a vertical bar with a curved top and a horizontal base. The 'E' is a vertical bar with a horizontal top bar and a horizontal middle bar. The 'L' at the bottom is a vertical bar with a horizontal top bar.