

Información general

Nombre de producto: EQUIPO DE FILTRACIÓN CON KITASATO, TAPÓN Y EMBUDO

Fecha Rev.: Jul 11, 2023

Detalles:

Compuesto por los siguientes productos:

- Embudo con placa filtrante, porosidad nº3 (Ref. 0100260563)
- Tapón de silicona con orificio (Ref. 0200259950)
- Matraz kitasato de 1000mL, con oliva de plástico (Ref. 0200913213) y vidrio (Ref. 0200913214)

Fotografía



Datos técnicos

Ficha de producto

Ref.	Oliva	Capacidad (ML)
0100259907	Plástico	1000
0100259908	Vidrio	1000



General Information

Product name: FILTRATION SET WITH FILTERING FLASK, STOPPER AND FUNNEL

Rev. Date: Jul 11, 2023

Details:

Includes the following products:

- Filter funnel with fritted disc, porosity nº3 (Ref. 0100260563)
- Silicone stopper with hole (Ref. 0200259950)
- Filter flask, 1000mL, with plastic (Ref. 0200913213) and glass (Ref. 0200913214) olive

Picture



Technical Information

Product datasheet

Ref.	Hose	Capacity (ML)
0100259907	Plastic	1000 ML
0100259908	Glass	1000 ML

A large, light green watermark logo for 'LORIMER' is oriented vertically on the right side of the page. The logo consists of a stylized 'L' at the top, followed by the letters 'O', 'R', 'I', 'M', 'E', 'R' in a bold, sans-serif font. Below the letters is a circular graphic element with a white outline, and at the bottom is a square graphic element with a white outline.