

## Información general

**Nombre de producto:** EQUIPO DE FILTRACIÓN CON ERLLENMEYER Y SOPORTE FILTRANTE CON OLIVA, UNIÓN ROSCADA TIPO SAV

**Fecha Rev.:** Jul 11, 2023

**Detalles:**

Compuesto por los siguientes productos:

- Cuerpo con placa porosa y oliva, macho (Ref. 0100259885)
- Embudo de filtración de 300mL (Ref. 0100260004)
- Pinza de sujeción (Ref. 0100260006)
- Matraz Erlenmeyer de 1000mL, con boca roscada de tipo SAV (Ref. 0200913314)
- Junta de nitrilo (Ref. 10JUNR0045)
- Tapón con orificio y rosca Rodaviss (Ref. 10TAOR0045)

## Fotografía



ERLENMEYER

Datos técnicos

Ref.	Descripción
0100259882	Equipo completo



## General Information

**Product name:** FILTRATION SET WITH ERLLENMEYER AND FILTER HOLDER WITH OLIVE, SAV TYPE THREAD JOINT

**Rev. Date:** Jul 11, 2023

**Details:**

Includes the following products:

- Body with fritted disc and olive, cone joint (Ref. 0100259885)
- Filter funnel, 300mL (Ref. 0100260004)
- Clamp (Ref. 0100260006)
- Flask, 1000mL, with SAV-type threaded neck (Ref. 0200913314)
- Nitrile seal (Ref. 10JUNR0045)
- Rodaviss screwcap, with aperture (Ref. 10TAOR0045)

## Picture



## Technical Information

Ref.	Description
0100259882	Complete set

L  
E  
E  
O  
L