

Información general

Nombre de producto: SOPORTE PARA PIPETAS

Fecha Rev.: Jul 11, 2023

Detalles:

Esterilizable, fabricado en PP.

El soporte para pipetas cuenta con dos placas giratorias, la inferior para proteger las puntas de las pipetas graduadas y aforadas, al mismo tiempo que retiene los restos de líquido que gotea.

Se suministra desmontado con instrucciones de montaje.

Fotografía



Datos técnicos

Ref.	Capacidad
10SOPT0010	94 pipetas

General Information

Product name: PIPETTE STAND

Rev. Date: Jul 11, 2023

Details:

Sterilizable, made from PP.

The pipette stand is comprised by two spinning plates. The lower one protects the graduated and on-mark pipette tips, while collecting the dripping liquids.

Supplied with assembly instructions.

Picture




Technical Information

Ref.	Capacity
10SOPT0010	94 pipettes

SOPORTE DE PIPETAS - PIPETTE STAND

PARTES DEL SOPORTE



- Base (1)
- Tubo largo (2)
- Escurreidor (3)
- Tubo corto (4)
- Rejilla (5)



ELEMENTS

- Base (1)
- Long tube (2)
- Wringer (3)
- Short tube (4)
- Rack (5)

Montaje/ Assembly



paso 3 / step 3



paso 4 / step 4

