

MPW-260

MPW-260R

MPW-260RH



MPW - 260

SPEED	4550	0	■
RDF	2486	0	■
TIME	09:02:00	00:02:00	

MPW MED. INSTRUMENTS

SHORT START STOP COVER BACK

PDF



FEATURES

effective ventilation system [MPW-260]

with cooling [MPW-260R]

with cooling and heating [MPW-260RH]

temperature +4°C ensured for maximum speed of each rotor [MPW-260R; MPW-260RH]

delayed start function

start when preselected temperature is reached [MPW-260R; MPW-260RH]

initial cooling without/with centrifuging, cooling after centrifuging, precise temperature stabilization [MPW-260R; MPW-260RH]

maintenance-free induction motor

modern software system

large graphic LCD [MPW-260RH] - blanking the screen after a period of inactivity

preselected and real values of speed, RCF, time shown on the display at the same time [MPW-260]

preselected and real values of speed, RCF, time and temperature shown on the display at the same time [MPW-260R; MPW-260RH]

100 user programs

10 acceleration/deceleration curves

programming of multi-section acceleration/deceleration curves

RPM/RCF adjustment

two modes of time calculation: from start key pressing or from reaching preselected speed

continuous operation mode - HOLD

SHORT time operation mode

manual centrifugal radius adjustment [with automatic RCF correction]

possibility of changing of parameters during centrifuging

automatic rotor identification

automatic lid opening

menu language setting [PL, EN, ES, IT, PT, DE, RU, SE, FR, CZ]

selected functions blocking, password protection

manual adjustment of density of samples >1,2g/cm³ [with automatic maximum speed correction]

stainless steel rotor chamber

SAFETY

unbalance sensor

lid locking during rotor running

start blocking at opened lid

emergency lid lock release

product conforming with the EN-61010-1 and EN-61010-2-020, EN-61010-2-101 standards

WYROBY MEDYCZNE DLA PROFESJONALISTÓW MEDICAL DEVICES FOR PROFESSIONAL USE ONLY



PN-EN ISO 9001
PN-EN ISO 13485



PRODUCENT WIRÓWEK LABORATORYNYCH
MANUFACTURER OF LABORATORY CENTRIFUGES

MPW MED. INSTRUMENTS

ul. Boremlowska 46 | 46 Boremlowska Street
04-347 Warszawa, Polska | 04-347 Warsaw, Polska

www.mpw.pl