

MPW-351e



PDF



v3



## MPW-351e

ⓘ 230V 50/60 Hz; 120V 50/60 Hz  
[opcja; [option](#); 100, 110, 127V 60Hz]

⚡ 360 W

🗑️ 800 ml

RPM max. 300 ÷ 4 500 RPM, krok/[step](#) 100 RPM

RCF max. 3 600 x g

🕒 1min. ÷ 99min. krok/[step](#) 1min.

✂️ 380 x 430 x 540mm [HxWxD]

📊 ~ 40,3/43,8 kg [230V/120V]

Cat. No. 10351e/2-56 [230V 50/60Hz]  
10351e/1-56 [120 50/60Hz]



### CECHY

bezobsługowy silnik indukcyjny

intuicyjny interfejs użytkownika

wyświetlacz LED

2 charakterystyki rozpędzania [FAST, SOFT]

4 charakterystyk hamowania [FAST, SOFT, LONG, RUN]

5 programów

praca w trybie SHORT

efektywny system wentylacji

szeroka gama wirników i wyposażenia

komora wirowania ze stali nierdzewnej

### BEZPIECZEŃSTWO

system wykrywania niewyważenia

blokada pokrywy podczas wirowania

blokada startu przy otwartej pokrywie

awaryjne otwieranie pokrywy

zabezpieczenie termiczne silnika

produkt zgodny z normami EN-61010-1  
i EN-61010-2-020, EN-61010-2-101

### FEATURES

maintenance-free induction motor

intuitive user interface

LED display

2 acceleration characteristics [FAST, SOFT]

4 deceleration characteristics [FAST, SOFT, LONG, RUN]

5 programs

SHORT mode operation

effective ventilation system

wide range of rotors and equipment

stainless steel rotor chamber

### SAFETY

unbalance detection system

lid locking during rotor running

start blocking at opened lid

emergency lid lock release

thermal motor protection

product conforming with the EN-61010-1  
and EN-61010-2-020, EN-61010-2-101 standards

### UWAGA

Zastrzegamy sobie prawo do zmian w materiałach reklamowych w dowolnym momencie w przypadku zmian w naszej ofercie online. MPW MED. INSTRUMENTS nie ponosi odpowiedzialności za ewentualne błędy drukarskie.

### PLEASE NOTE

We reserve the right to change information included in our promotional materials at any time in the case the online offer changes. MPW MED. INSTRUMENTS is not responsible for printing mistakes.

# WYROBY MEDYCZNE DLA PROFESJONALISTÓW MEDICAL DEVICES FOR PROFESSIONAL USE ONLY



PRODUCENT WIRÓWEK LABORATORYNYCH  
MANUFACTURER OF LABORATORY CENTRIFUGES

**MPW MED. INSTRUMENTS**

ul. Boremlowska 46 | 46 Boremlowska Street  
04-347 Warszawa, Polska | 04-347 Warsaw, Polska

[www.mpw.pl](http://www.mpw.pl)