


# NEUATION



DESIGNED TO  
DELIVER

**iStir OH Jr.**  
(Analog Overhead Stirrer)

## iStir OH Jr. (Analog Overhead Stirrer)

Overhead stirrers are a crucial part of laboratory setups and the work that they do makes the user's life a lot easier. But, in the race to attract customers most manufacturers skimp on the really important stuff and focus on the superficial to sell a lot of devices. In the end, it's the customer who ends up suffering because most of these devices would suffer from some kind of performance or maintenance issues.

Here at Neuation, we start from the basics to create a product that not only performs well, but also creates a pleasant experience for the user. With a complete department focusing on R&D, the designers and product developers at neuation are armed with crucial research about the customer's requirements and needs when it comes to overhead stirrers. The result as you can see is the OH Jr. series of stirrers which combine a compact form factor with exceptional performance for enhanced convenience. Various built-in safety features ensure that no matter what, your safety is never compromised.



## Product Features



Brushless DC motor delivers ultra quiet performance and a maintenance-free long life



Variable stirring speed settings upto 2200 RPM for powerful performance



Small footprint saves valuable bench space

IP42 class compliant rating for unmatched durability



Built-in overload safety feature to ensure that mindless mistakes do not compromise device & your safety



Powerful motor delivers torque upto 10 Ncm to stir solutions with viscosity up to 1000 mPas & 10 L volume



## Specifications

iStir OH Jr.	
Specifications	Details
Motor Type	Brushless DC
Capacity	10 Liters
Speed Range	60 - 2200 RPM
Torque	10 Ncm
Vascocity	1000 mPas
Timer	1 min to 99 hr 59 mins & infinity mode
Protection Class	IP 42
Dimension (Head)	80 x 42 x 166mm
Weight	1.3 kgs

## Ordering Information

iStir OH Jr.	
Product Name	Cat #
iStir OH Jr.	OHS1-1100-EU/US/UK/AU

Accessories		
Product Name	Cat #	
Boss Head	OHS1-BH	
Impellers	Length (mm)	
OHS1-JX	300	Jr. impeller

# NEUATION

www.neuation.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.